

# Introducing the **enovate**<sup>™</sup> medical carts



Enovate Medical Laptop Cart

Enovate Medical LCD Cart

Clean. Simple. Smart.  
**Mobility just got easier.**

MADE IN USA

MADE IN USA



## » SMART TECHNOLOGY

A smart dashboard, located at the caregiver's fingertips, provides immediate access to the battery alert, power on/off, LED fuel gauge, keyboard light and computer reboot.\*

## » ERGONOMICS

Designed to be comfortable & ergonomically correct whether you are pushing, pulling, sitting, or standing. Enovate's keyboard, monitor, and work surfaces can be adjusted independently for maximum flexibility and clinician comfort.



## » POWER OPTIONS

Multiple options for power, from NiMH to a hot-swappable Lithium-Ion. Simple maintenance, provides quick and easy accessibility to the power system through a drawer-style compartment.

### Type

- Integrated
- Hot-Swap
- Non-powered

### Current

- AC
- DC

### Features

- Certified to UL 60601-1
- Output: Pure Sine Wave
- Gentle shut-down at critical low power levels
- Automatic Transfer Switch enables full use of computer without affecting charge time.

### Chemistry


- SLA
- NiMH
- Lithium-Ion

### Capacity

- 20 Amp Hour
- 40 Amp Hour
- 50 Amp Hour
- 200 Watt
- 280 Watt
- 360 Watt

### WARRANTY

Standard Warranty: 4-year on cart; 2-year on electrical; 6 month limited warranty on battery.



## » INFECTION CONTROL

All contact surfaces are finished or made with anti-microbial materials delivering smooth and comfortable surfaces.

### STANDARD SPECIFICATIONS

- Antimicrobial Contact Surfaces
- Smart Dashboard
  - » Battery Alert
  - » Power On/Off
  - » LED Fuel Gauge
  - » Remote PC Reboot Button\*
- Keyboard Light with Auto Shut-off
- Mouse Bracket
- Left/Right Mouse Tray
- Comfort Wrist Rest
- Comfort Grip Front Handle
- Comfort Grip Rear Handle
- Comfort Grip Height Adjustment Lever
- 17" x 17" Work Surface
- Cable Storage
- 3 port USB Hub
- Internal Cable Management
- Ethernet Ready
- 4 casters (2 locking)
- Hospital Grade Spiral Cord
- Keyed Lock Security
- Power Management Software
- LCD Height Adjustment: simple landscape to portrait adjustment
- Integrated Power Conditioning System
- Internal CPU Storage: 16" w x 16" d x 2.75" h
- Negative-tilt Keyboard
- Weight:
  - » Powered: Starting as low as 100 lbs!
  - » Non-Powered: 65 lbs

\* Enovate Medical Cart works with selected computer devices.

\*\* Enovate Medical Cart power system has some components that are imported.

enovate™

1250 Woodward Heights  
Ferndale, Mi 48220-1427  
p. (877)258-8030 f. (847)890-6131

[www.enovateusa.com](http://www.enovateusa.com)  
[salesrequest@enovateit.com](mailto:salesrequest@enovateit.com)

©2009 EnovateIT LLC. All Rights Reserved

enovate™