

Enovate Medication Carts

Less trips to the medication room means higher efficiency and fewer errors for clinicians to worry about. We've delivered essential features in a compact form ensuring medication and supplies stay secure and organized at the point-of-care.



Enovate Laptop Medication Cart

Enovate LCD Medication Cart

Simple. Smart. Lightweight.
Med-dispensing just got easier.

DRAWER BINS

**Small Bin**

Interior Dimension
3.125" W x 2.125" H x 10.5" D
(825 mm x 571 mm x 265 mm)

**Medium Bin**

Interior Dimension
7.25" W x 2.125" H x 10.5" D
(184 mm x 571 mm x 264 mm)

Robot-Ready Option Available

**Large Bin**

Inside Dimensions
14" w x 2.5" h x 10.25" d
(260 mm x 64 mm x 274 mm)

Optional Five Dividers, Six Spaces (Left to Right) 2.22" (56 mm) Spacing Available

**Small 2-Tier Deep Bin**

Interior Dimension
3.125" W x 5.25" H x 10.5" D
(825 mm x 571 mm x 265 mm)

**Medium 4-Tier Deep Bin**

Envelope-Ready Envelope Friendly
Inside Dimensions
6.75" w x 9.5" h x 10.5" d
(184 mm x 241 mm x 264 mm)

Optional Four Dividers, Five 1.25" (318 mm) Spaces for Envelopes Available

**Large Deep Bin**

Inside Dimensions
14" w x 4.22" h x 10.77" d
(356 mm x 107 mm x 274 mm)

Optional Five Dividers, Six Spaces (Left to Right) 2.22" (56 mm) Spacing Available

**Large 4-Tier Deep Bin**

Inside Dimensions
14" w x 10" h x 10.77" d
(356 mm x 254 mm x 274 mm)

Optional Five Dividers, Six Spaces (Front to Back) 1.79" (46 mm) Spacing Available

Designed with superior capacity and modularity — choose from multiple drawer configurations and bin sizes to fit your specific requirements.

POWER OPTIONS†

Multiple options for power. Simple maintenance. Quick and easy accessibility to the power system through a drawer-style compartment.

Chemistry

- EON Lithium Phosphate
- Sealed Lead Acid (SLA)

Currents

- AC
- DC

Capacity

- 20 Amp Hour
- 40 Amp Hour

Features

- Certified to UL 60601-1
- Gentle shut-down at critical low power levels.
- Automatic Transfer Switch enables full use of computer without affecting charge time.

Warranty

Standard Warranty: four years on structural; two years on electrical; six month limited warranty on SLA batter; three years limited warranty on EON Lithium Phosphate battery.

AVAILABLE NOW!



SMART TECHNOLOGY

An intelligent control panel is located at the caregiver's fingertips providing immediate access to the cart's battery alert, battery power button, battery LED fuel gauge, a work light and a laptop reboot button for select CPUs/Laptops (Please contact your local representative to find out what CPUs/laptops are compatible).

The Enovate Medication Carts offer an array of safety features, including anti-microbial medication drawers and supply bins that are protected by a keyless electronic PIN code and a programmable auto lock.

Software

Remote pin code management, power monitoring and battery management software included.

Wireless

Wireless-friendly, laptop/CPU tray was designed to optimize wireless propagation.

Security

All technology, supplies and medication are securely stored to help prevent theft and misuse through keyed and electronic locks.

SPECIFICATIONS

All contact surfaces are finished or made with anti-microbial materials and are designed to be smooth, comfortable and easy to clean —to help control infection and prevent dangerous bacteria from spreading.

Cart Specifications

- Made in USA
- Antimicrobial Contact Surfaces
- Smart Control Panel
 - » Battery Alert
 - » Power On/Off
 - » LED Fuel Gauge
 - » Remote PC Reboot Button‡
- Drawer Keypad: Electronic PIN codes provide keyless entry to medication and supplies.
- Work Light with Auto Shut-off
- Mouse Bracket/Holder
- Mousing Surface: Left/Right Mouse Tray
- Work Surface that helps control spills 17" x 17" (432 mm x 432 mm)

- Cable Storage
- 3 port USB Hub
- Internal Cable Management
- Ethernet Ready
- 4 Durable Precision Casters (2 Locking)
- Hospital Grade Spiral Power Cord
- Cord Hanger: used to store the power cord within reach
- LCD Height Adjustment: Simple landscape to portrait adjustment
- Integrated Power Conditioning System
- Internal CPU Storage: 15.875" w x 11.875" d x 2.875" h (403 mm x 302 mm x 73 mm)

Drawer Configurations

- Multiple Drawer Configurations Available: Drawers are available in configurations from 1-16 bins.
- Infection Control: Anti-microbial powder coating & plastics
- Bin/Patient Labeling: Identification on drawer bin face
- Removable Tilt Locking Bins supported
- Drawer Keypad: Electronic PIN codes provide keyless entry to medications and supplies
- Includes an Emergency Keyed Lock and Override Switch (One key included)
- Customizable auto-lock timer for drawers

- Bin Dividers
- Multiple User Pin Codes
- Upgradeable for the Future: Future smart drawers that open individually per patient.

Standard Accessories on Medication Dispensing Cart Bundles

- Plastic Overlay Document Protector
- Utility Bin/Waste Basket
- Scanner Bracket
- Sharps Bracket



enovate™

7820 N. Lilley Road
Canton, MI 48187
p. (877)258-8030 f. (847)890-6131

www.enovateusa.com
salesrequest@enovateit.com

† Enovate Medical Cart works with selected computer devices.

‡ Enovate Medical Cart power system has some components that are imported.